

Product Datasheet

ATP2A1 peptide (orb216947)

Description	The peptide is used to block anti-SERCA1 antibody reactivity.
Reactivity	Human, Mouse, Rabbit, Rat
Tested Applications	Blocking
Target	ATP2A1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213591.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q8R429, Q64578, O14983
Entrez	116601, 11937, 487
Expiration Date	12 months from date of receipt.