

Product Datasheet

ATP6AP1 peptide (orb216950)

Description	The peptide is used to block anti-ATP6AP1 antibody reactivity.
Reactivity	Human, Mouse, Porcine, Rat
Tested Applications	Blocking
Target	ATP6AP1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213594.
Purity	> 85%
Source	Synthetic
Uniprot ID	O54715 , Q9R1Q9 , Q15904
Entrez	83615 , 54411 , 537
Expiration Date	12 months from date of receipt.