

Product Datasheet

BAX peptide (orb216955)

Description	BAX peptide
Reactivity	Human, Mouse, Porcine, Rat
Tested Applications	Blocking
Target	BAX
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213599.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q63690, Q07812
Entrez	581
Expiration Date	12 months from date of receipt.