

Product Datasheet

BCL2 peptide (orb216958)

Description	BCL2 peptide
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	BCL2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213602.
Purity	> 85%
Source	Synthetic
Uniprot ID	P10415, P49950, P10417
Entrez	12043, 596, 24224
Expiration Date	12 months from date of receipt.