

Product Datasheet

BMP7 peptide (orb216964)

Description	BMP7 peptide
Reactivity	Human, Mouse, Porcine, Primate, Rabbit, Rat
Tested Applications	Blocking
Target	BMP7
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213608.
Purity	> 85%
Source	Synthetic
Uniprot ID	P23359, P18075
Entrez	655, 12162
Expiration Date	12 months from date of receipt.