

Product Datasheet

BRCA1 peptide (orb216969)

Description	BRCA1 peptide
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	BRCA1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213613.
Purity	> 85%
Source	Synthetic
Uniprot ID	P48754, O54952, P38398
Entrez	497672, 672, 12189
Expiration Date	12 months from date of receipt.