

Product Datasheet

ZFP36L1 peptide (orb216972)

Description	ZFP36L1 peptide
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	ZFP36L1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213616.
Purity	> 85%
Source	Synthetic
Uniprot ID	P23950, P17431, Q07352
Entrez	29344, 677, 12192
Expiration Date	12 months from date of receipt.