

Product Datasheet

CA1 peptide (orb216982)

Description	The peptide is used to block anti-Carbonic Anhydrase 1 antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	CA1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213626.
Purity	> 85%
Source	Synthetic
Uniprot ID	B0BNN3, P13634, P00915
Entrez	12346, 759, 310218
Expiration Date	12 months from date of receipt.