

Product Datasheet

CA2 peptide (orb216983)

Description	The peptide is used to block anti-Carbonic Anhydrase 2 antibody reactivity.
Reactivity	Human, Primate, Rat
Tested Applications	Blocking
Target	CA2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213627.
Purity	> 85%
Source	Synthetic
Uniprot ID	P27139, P00918
Entrez	760, 54231
Expiration Date	12 months from date of receipt.