

Product Datasheet

CACNG1 peptide (orb216988)

Description	The peptide is used to block anti-CACNG1 antibody reactivity.
Reactivity	Human, Mouse, Porcine, Rabbit, Rat
Tested Applications	Blocking
Target	CACNG1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213632.
Purity	> 85%
Source	Synthetic
Uniprot ID	P97707, Q06432, O70578
Entrez	786, 12299, 29658
Expiration Date	12 months from date of receipt.