

Product Datasheet

CAMK2B peptide (orb216990)

Description	The peptide is used to block anti-CaMK2 beta antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	CAMK2B
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213634.
Purity	> 85%
Source	Synthetic
Uniprot ID	P28652, Q13554, P08413
Entrez	816, 12323, 24245
Expiration Date	12 months from date of receipt.