

Product Datasheet

Caspase 1 peptide (orb216995)

Description	Caspase 1 peptide
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	CASP1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213639.
Purity	> 85%
Source	Synthetic
Uniprot ID	P29466, P29452, P43527
Entrez	12362, 834
Expiration Date	12 months from date of receipt.