

Product Datasheet

CASP2 peptide (orb216998)

Description	CASP2 peptide
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	CASP2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213642.
Purity	> 85%
Source	Synthetic
Uniprot ID	P29594, P55215, P42575
Entrez	835, 12366, 64314
Expiration Date	12 months from date of receipt.