

Product Datasheet

Caspase 3 peptide (orb217002)

Description	Caspase 3 peptide
Reactivity	Human, Mouse, Porcine, Primate, Rat
Tested Applications	Blocking
Target	CASP3
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213646.
Purity	> 85%
Source	Synthetic
Uniprot ID	P70677, P42574
Entrez	836, 12367
Expiration Date	12 months from date of receipt.