

## Product Datasheet

### CASP6 peptide (orb217011)

|                            |                                                          |
|----------------------------|----------------------------------------------------------|
| <b>Description</b>         | CASP6 peptide                                            |
| <b>Reactivity</b>          | Human, Mouse, Rat                                        |
| <b>Tested Applications</b> | Blocking                                                 |
| <b>Target</b>              | CASP6                                                    |
| <b>Form/Appearance</b>     | Lyophilized powder                                       |
| <b>Storage</b>             | Shipped at 4°C. Store at -20°C for one year.             |
| <b>Note</b>                | For research use only                                    |
| <b>Application notes</b>   | This peptide was used as immunogen to produce orb213655. |
| <b>Purity</b>              | > 85%                                                    |
| <b>Source</b>              | Synthetic                                                |
| <b>Uniprot ID</b>          | <b>O35397</b> , <b>P55212</b> , <b>O08738</b>            |
| <b>Entrez</b>              | <b>12368</b> , <b>839</b>                                |
| <b>Expiration Date</b>     | 12 months from date of receipt.                          |