

Product Datasheet

CASP7 peptide (orb217014)

Description	CASP7 peptide
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	CASP7
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213658.
Purity	> 85%
Source	Synthetic
Uniprot ID	P97864, P55210
Entrez	12369, 840
Expiration Date	12 months from date of receipt.