

Product Datasheet

Caspase 9 peptide (orb217018)

Description	Caspase 9 peptide
Reactivity	Human, Primate
Tested Applications	Blocking
Target	CASP9
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213662.
Purity	> 85%
Source	Synthetic
Uniprot ID	P55211
Entrez	842
Expiration Date	12 months from date of receipt.