

Product Datasheet

Caveolin 1 peptide (orb217023)

Description	Caveolin 1 peptide
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	CAV1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213667.
Purity	> 85%
Source	Synthetic
Uniprot ID	P41350, P49817, Q03135
Entrez	857, 12389
Expiration Date	12 months from date of receipt.