

Product Datasheet

CBL peptide (orb217028)

Description	The peptide is used to block anti-c-CBL
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	CBL
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213672.
Purity	> 85%
Source	Synthetic
Uniprot ID	P22681, P22682
Entrez	12402, 867
Expiration Date	12 months from date of receipt.