

Product Datasheet

CCNF peptide (orb217035)

Description	The peptide is used to block anti-Cyclin F antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	CCNF
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213679.
Purity	> 85%
Source	Synthetic
Uniprot ID	P41002, P51944, Q8K4F8
Entrez	117524, 899, 12449
Expiration Date	12 months from date of receipt.