

Product Datasheet

41889 peptide (orb217054)

Description	The peptide is used to block anti-Septin 7 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	SEPT7
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213698.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q16181 , O55131 , Q9WVC0
Entrez	64551 , 989 , 235072
Expiration Date	12 months from date of receipt.