

Product Datasheet

N Cadherin peptide (orb217063)

Description	N Cadherin peptide
Reactivity	Human, Mouse, Rat, Zebrafish
Tested Applications	Blocking
Target	CDH2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213707.
Purity	> 85%
Source	Synthetic
Uniprot ID	P15116, P19022, Q9Z1Y3
Entrez	83501, 12558, 1000
Expiration Date	12 months from date of receipt.