

Product Datasheet

CEBPA peptide (orb217078)

Description	The peptide is used to block anti-C/EBP alpha
Reactivity	Human, Mouse, Rat, Zebrafish
Tested Applications	Blocking
Target	CEBPA
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213722.
Purity	> 85%
Source	Synthetic
Uniprot ID	P05554, P53566, P49715
Entrez	24252, 12606, 1050
Expiration Date	12 months from date of receipt.