

## Product Datasheet

### CEBPB peptide (orb217079)

|                            |                                                                   |
|----------------------------|-------------------------------------------------------------------|
| <b>Description</b>         | The peptide is used to block anti-C/EBP beta antibody reactivity. |
| <b>Reactivity</b>          | Human, Mouse, Porcine, Rat                                        |
| <b>Tested Applications</b> | Blocking                                                          |
| <b>Target</b>              | CEBPB                                                             |
| <b>Form/Appearance</b>     | Lyophilized powder                                                |
| <b>Storage</b>             | Shipped at 4°C. Store at -20°C for one year.                      |
| <b>Note</b>                | For research use only                                             |
| <b>Application notes</b>   | This peptide was used as immunogen to produce orb213723.          |
| <b>Purity</b>              | > 85%                                                             |
| <b>Source</b>              | Synthetic                                                         |
| <b>Uniprot ID</b>          | <b>P17676, P21272, P28033</b>                                     |
| <b>Entrez</b>              | <b>24253, 1051, 12608</b>                                         |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                   |