

Product Datasheet

CHGA peptide (orb217092)

Description	The peptide is used to block anti-Chromogranin A antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	CHGA
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213736.
Purity	> 85%
Source	Synthetic
Uniprot ID	P10645, P10354, P26339
Entrez	12652, 1113
Expiration Date	12 months from date of receipt.