

Product Datasheet

CHRNA7 peptide (orb217096)

Description	The peptide is used to block anti-CHRNA7 antibody reactivity.
Reactivity	Human, Mouse, Porcine, Primate, Rat, Zebrafish
Tested Applications	Blocking
Target	CHRNA7
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213740.
Purity	> 85%
Source	Synthetic
Uniprot ID	P49582 , Q05941 , Q494W8 , P36544
Entrez	25302 , 11441 , 1139
Expiration Date	12 months from date of receipt.