

Product Datasheet

CHUK peptide (orb217097)

Description	The peptide is used to block anti-IKK alpha
Reactivity	Human, Mouse, Porcine, Rat
Tested Applications	Blocking
Target	CHUK
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213741.
Purity	> 85%
Source	Synthetic
Uniprot ID	O15111, Q60680
Entrez	12675, 1147
Expiration Date	12 months from date of receipt.