

Product Datasheet

CKM peptide (orb217101)

Description	The peptide is used to block anti-Creatine Kinase M antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	CKM
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213745.
Purity	> 85%
Source	Synthetic
Uniprot ID	P07310, P00564, P06732
Entrez	24265, 12715, 1158
Expiration Date	12 months from date of receipt.