

Product Datasheet

AP2M1 peptide (orb217102)

Description	The peptide is used to block anti-AP2M1 antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat, Zebrafish
Tested Applications	Blocking
Target	AP2M1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213746.
Purity	> 85%
Source	Synthetic
Uniprot ID	P84092, Q96CW1, P84091
Entrez	116563, 11773, 1173
Expiration Date	12 months from date of receipt.