

Product Datasheet

Collagen I peptide (orb217113)

Description	Collagen I peptide
Reactivity	Human, Mouse, Porcine, Rat, Sheep
Tested Applications	Blocking
Target	COL1A2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213757.
Purity	> 85%
Source	Synthetic
Uniprot ID	P02466, Q01149, P08123
Entrez	1278, 84352, 12843
Expiration Date	12 months from date of receipt.