

Product Datasheet

Collagen I peptide (orb217114)

Description	Collagen I peptide
Reactivity	Human, Rabbit
Tested Applications	Blocking
Target	COL1A2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213758.
Purity	> 85%
Source	Synthetic
Uniprot ID	P08123
Entrez	1278
Expiration Date	12 months from date of receipt.