

Product Datasheet

Collagen III peptide (orb217115)

Description	Collagen III peptide
Reactivity	Human, Mouse, Porcine, Rabbit, Rat
Tested Applications	Blocking
Target	COL3A1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213759.
Purity	> 85%
Source	Synthetic
Uniprot ID	P02461, P13941, P08121
Entrez	84032, 12825, 1281
Expiration Date	12 months from date of receipt.