

Product Datasheet

Collagen VIII peptide (orb217120)

Description	Collagen VIII peptide
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	COL8A1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213764.
Purity	> 85%
Source	Synthetic
Uniprot ID	P27658, Q00780
Entrez	12837, 1295
Expiration Date	12 months from date of receipt.