

Product Datasheet

COL11A1 peptide (orb217121)

Description	The peptide is used to block anti-Collagen 11 alpha 1 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	COL11A1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213765.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q61245 , P12107 , P20909
Entrez	25654 , 12814 , 1301
Expiration Date	12 months from date of receipt.