

Product Datasheet

COMT peptide (orb217122)

Description	The peptide is used to block anti-COMT antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	COMT
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213766.
Purity	> 85%
Source	Synthetic
Uniprot ID	P21964, P22734, O88587
Entrez	24267, 12846, 1312
Expiration Date	12 months from date of receipt.