

Product Datasheet

CPA1 peptide (orb217123)

Description	The peptide is used to block anti-Carboxypeptidase A1 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	CPA1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213767.
Purity	> 85%
Source	Synthetic
Uniprot ID	P00731, P15085, Q7TPZ8
Entrez	109697, 1357, 24269
Expiration Date	12 months from date of receipt.