

Product Datasheet

CPE peptide (orb217124)

Description	The peptide is used to block anti-Carboxypeptidase E antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	CPE
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213768.
Purity	> 85%
Source	Synthetic
Uniprot ID	P15087, P16870, Q00493
Entrez	1363, 25669, 12876
Expiration Date	12 months from date of receipt.