

Product Datasheet

CREB1 peptide (orb217133)

Description	CREB1 peptide
Reactivity	Human, Mouse, Porcine, Rat, Sheep, Zebrafish
Tested Applications	Blocking
Target	CREB1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213777.
Purity	> 85%
Source	Synthetic
Uniprot ID	P15337, P16220, Q01147
Entrez	81646, 12912, 1385
Expiration Date	12 months from date of receipt.