

Product Datasheet

ATF2 peptide (orb217135)

Description	ATF2 peptide
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	ATF2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213779.
Purity	> 85%
Source	Synthetic
Uniprot ID	P15336, Q00969, P16951
Entrez	11909, 81647, 102641666, 1386
Expiration Date	12 months from date of receipt.