

Product Datasheet

CREBBP peptide (orb217136)

Description	CREBBP peptide
Reactivity	Human, Mouse, Porcine, Rat
Tested Applications	Blocking
Target	CREBBP
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213780.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q6JHU9 , Q92793 , P45481
Entrez	1387
Expiration Date	12 months from date of receipt.