

Product Datasheet

CRP peptide (orb217138)

Description	The peptide is used to block anti-C-Reactive Protein antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	CRP
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213782.
Purity	> 85%
Source	Synthetic
Uniprot ID	P48199, P14847, P02741
Entrez	25419, 1401, 12944
Expiration Date	12 months from date of receipt.