

Product Datasheet

CSE1L peptide (orb217140)

Description	The peptide is used to block anti-XPO2 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	CSE1L
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213784.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9ERK4, P55060
Entrez	110750, 1434
Expiration Date	12 months from date of receipt.