

Product Datasheet

CSF2RB peptide (orb217141)

Description	The peptide is used to block anti-CD131 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	CSF2RB
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213785.
Purity	> 85%
Source	Synthetic
Uniprot ID	P26955, P32927
Entrez	12983, 1439
Expiration Date	12 months from date of receipt.