

Product Datasheet

CSNK1G2 peptide (orb217146)

Description	The peptide is used to block anti-CK1 gamma 2 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	CSNK1G2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213790.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q8BVP5 , Q62762 , P78368
Entrez	103236 , 1455 , 65278
Expiration Date	12 months from date of receipt.