

Product Datasheet

CSNK1G3 peptide (orb217147)

Description	The peptide is used to block anti-CK1 gamma 3 antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	CSNK1G3
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213791.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q62763, Q9Y6M4, Q8C4X2
Entrez	1456, 64823, 70425
Expiration Date	12 months from date of receipt.