

Product Datasheet

CSRP1 peptide (orb217151)

Description	The peptide is used to block anti-CSRP1 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	CSRP1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213795.
Purity	> 85%
Source	Synthetic
Uniprot ID	P97315, P21291, P47875
Entrez	13007, 29276, 1465
Expiration Date	12 months from date of receipt.