

Product Datasheet

CXADR peptide (orb217167)

Description	The peptide is used to block anti-CXADR antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	CXADR
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213811.
Purity	> 85%
Source	Synthetic
Uniprot ID	P97792, Q9R066, P78310
Entrez	89843, 1525, 13052
Expiration Date	12 months from date of receipt.