

Product Datasheet

CYP2A6 peptide (orb217175)

Description	The peptide is used to block anti-Cytochrome P450 2A6 antibody reactivity.
Reactivity	Human, Primate
Tested Applications	Blocking
Target	CYP2A6; CYP2A7; CYP2A13
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213819.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q16696 , P11509 , P20853
Entrez	1553
Expiration Date	12 months from date of receipt.