

Product Datasheet

CYP2B6 peptide (orb217176)

Description	The peptide is used to block anti-Cytochrome P450 2B6 antibody reactivity.
Reactivity	Human
Tested Applications	Blocking
Target	CYP2B6
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213820.
Purity	> 85%
Source	Synthetic
Uniprot ID	P20813
Entrez	1555
Expiration Date	12 months from date of receipt.