

Product Datasheet

CYP2D6 peptide (orb217180)

Description	The peptide is used to block anti-Cytochrome P450 2D6 antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	CYP2D6
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213824.
Purity	> 85%
Source	Synthetic
Uniprot ID	P10635
Entrez	1565
Expiration Date	12 months from date of receipt.